

# Accessories Catalog



Discover your  
perfect **fit** from **CAIRE.**

# Table of Contents



## Item Category Page

<b>Oxygen Concentrator Accessories .....</b>	<b>3</b>
<b>Concentrator Batteries, Power Supplies, Backpacks, Bags, Carts, Belts .....</b>	<b>4</b>
AirSep External Power Cartridge & AirBelt .....	4
AirSep Focus.....	5
AirSep FreeStyle .....	7
AirSep FreeStyle 5 .....	8
AirSep LifeStyle .....	10
SeQual Eclipse.....	11
SeQual eQuinox .....	12
<b>Concentrator Therapeutic Accessories .....</b>	<b>14</b>
<b>Concentrator Testing and Service Equipment .....</b>	<b>15</b>
AirSep Focus, FreeStyle, FreeStyle 5 .....	15
AirSep NewLife and VisionAire Families .....	15
SeQual Eclipse.....	17
SeQual eQuinox .....	17
<b>Liquid Oxygen Accessories.....</b>	<b>18</b>
<b>Liquid Oxygen Portable Backpacks .....</b>	<b>19</b>
Companion .....	19
HELiOS .....	19
Hi Flow Stroller, Spirit, Sprint.....	19
<b>Liquid Oxygen Portable Hip Packs .....</b>	<b>20</b>
Companion .....	20
HELiOS .....	20
Spirit .....	21
<b>Liquid Oxygen Portable Universal Cart .....</b>	<b>22</b>
<b>Liquid Oxygen Reservoir Rollerbase Assemblies .....</b>	<b>23</b>
<b>Nasal Cannulas .....</b>	<b>24</b>
<b>Transfilling Accessories.....</b>	<b>25</b>
<b>Testing Accessories .....</b>	<b>29</b>
<b>General Accessories.....</b>	<b>34</b>
<b>Flow Upgrades and Conversions .....</b>	<b>35</b>
<b>CAIRE Generation 4 Meter Upgrades .....</b>	<b>38</b>
<b>Spirit Upgrade Kits .....</b>	<b>39</b>

NOTE: This catalog is not an exhaustive list of all of our available accessories. If items not included in this catalog are needed, please contact the Technical Support Team in your region for further assistance.

**Technical Support USA**  
800.482.2473  
techservice.usa@chartindustries.com

**Technical Support Europe**  
+44 (0) 1189 367060  
techservice.europe@chartindustries.com

# Oxygen Concentrator Accessories



Discover your  
perfect**fit** from CAIRE.







# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

AirSep External Power Cartridge & AirBelt™		
Item Image	Item Number	Description
	<i>Country-Specific Power Cords Available:</i> BT026-1 (USA) BT026-2 (Europe) BT026-3 (Australia) BT026-4 (UK) BT026-5 (Brazil) BT026-6 (Swiss)	<b>AirSep External Power Cartridge Kit:</b> Includes one external power cartridge (BT024-1), one charger (PW029), and power supply for the charger. For use with FreeStyle and FreeStyle 5. Also compatible with the Focus using the Focus Cable Adapter, PN: CD035-1. Please see below.
	BT024-1	<b>AirSep External Power Cartridge:</b> Extended battery power source for AirSep portable oxygen concentrators. Light-weight with built-in belt clip. Compatible with all AirSep portable oxygen concentrators including: Focus, FreeStyle, and FreeStyle 5.
	<i>Country-Specific Power Cords Available:</i> PW029-1 (USA) PW029-2 (Europe) PW029-3 (Australia) PW029-4 (UK) PW029-5 (Brazil) PW029-6 (Swiss)	<b>AirSep External Power Cartridge Charger:</b> Desktop charger for AirSep External Power Cartridges (BT024-1). Ability to charge 2 power cartridges simultaneously. Includes power supply for charging dock.
	<i>Country-Specific Power Cords Available:</i> BT017-1 (USA) BT017-2 (Europe) BT017-3 (Australia) BT017-4 (UK)	<b>AirBelt Battery:</b> Extended battery power source for AirSep portable oxygen concentrators designed to be worn around the waist. Compatible with all AirSep portable oxygen concentrators including: Focus, FreeStyle, and FreeStyle 5. (Includes power supply for recharging.)
	CD035-1	<b>Focus Cable Adapter for External Power Cartridge or AirBelt:</b> Cable designed specifically for connection of the Focus to the External Power Cartridge or AirBelt. Please note: FreeStyle and FreeStyle 5 portable oxygen concentrators do not require this item to connect to the External Power Cartridge or AirBelt.
	<i>Country-Specific Power Cords Available:</i> PW008-1S (USA) PW008-2S (Europe) PW008-3S (Australia) PW008-4S (UK)	<b>AC Power Supply for AirBelt:</b> Power Supply designed to recharge the AirBelt in a standard wall outlet. Includes power cord.








# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

AirSep Focus™		
Item Image	Item Number	Description
	BT023-1	<b>Focus Battery:</b> Additional battery for Focus portable oxygen concentrator.
	CD023-2	<b>AC Power Cord:</b> Detachable hospital-grade US power supply that connects to the AC Power Supply of the LifeStyle, FreeStyle 3, FreeStyle 5, or Focus. May be removed and interchanged with country-specific cords for travel. Contact Customer Service for country-specific cord part numbers.
	CD032-1	<b>Focus Battery Operation Coiled Cable (on/off switch):</b> Coil cord (with on/off switch) for the connection of the battery to the Focus to operate unit.
	CD032-2	<b>Focus Battery Charging Coiled Cable (no on/off switch):</b> Coil cord designed for the connection of the Focus battery to the Universal AC/DC Power Supply for recharging.
	CD034-1	<b>DC Auto Adaptor Cable for FreeStyle 3 and Focus:</b> DC cable to connect to the Universal Power Supply of the FreeStyle and Focus. Allows for charging and operation of the unit in an automobile.
	MI017-1	<b>FreeStyle &amp; Focus Carry-All Bag Strap:</b> Designed to allow the accessory bag to be carried over the shoulder. Attaches to MI320-1.
	MI320-1	<b>FreeStyle &amp; Focus Carry-All Bag:</b> Accessory bag designed specifically for easy transport of the FreeStyle and Focus portable oxygen concentrators and their accessories.

# Accessories

## Concentrator Batteries, Power Supplies, Backpacks, Bags, Carts, Belts

AirSep Focus <sup>™</sup>		
Item Image	Item Number	Description
Please see the figure below.	MI332-1	<b>Focus Bag Set:</b> Complete kit for the ultimate traveler! The Focus Bag Set includes the Battery Bag (MI333-1), Carrying Bag (MI345-1), Bag with Battery Pockets (MI379-1), and the Shoulder Strap (MI380-1) for extra comfort.
		
	MI333-1	<b>Focus Battery Case:</b> Battery case with belt clip designed to hold Focus battery.
	MI345-1	<b>Focus Carrying Bag with Belt Clip:</b> Carrying bag for Focus that contains a built-in belt clip. Also allows for attachment of the shoulder strap (not included).
	MI379-1	<b>Focus Carrying Bag with Battery Pockets:</b> Cloth bag designed specifically for the Focus that includes battery compartments for easy transportation of Focus and two batteries.
	MI380-1	<b>Focus Shoulder Strap:</b> Designed to allow the Focus to be carried over the shoulder. Can be attached to either MI379-1 or MI345-1 bags pictured above.
	PW019-1	<b>Universal Power Supply (Brick Only) For Focus:</b> Does not include any cords or cables.
	Country-Specific Power Cords Available: PW023-1 (USA) PW023-2 (Europe) PW023-3 (Australia) PW023-4 (UK)	<b>Universal Power Supply for Focus:</b> A complete kit which includes everything needed to operate and recharge the Focus on either AC or DC power. Includes items CD034-1 and CD023-2. Battery charging cable (CD032-2) is not included.



# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

AirSep FreeStyle™		
Item Image	Item Number	Description
	20678526	<b>Wheelchair Straps:</b> Straps designed for easy attachment of concentrators to wheelchairs. Can be used with the Eclipse, eQuinox, FreeStyle, and FreeStyle 5.
	BT018-1S	<b>FreeStyle 3 Battery:</b> Internal battery for the FreeStyle 3 portable oxygen concentrator.
	CD023-2	<b>AC Power Cord:</b> Detachable hospital-grade US power supply that connects to the AC Power Supply of the LifeStyle, FreeStyle 3, FreeStyle 5, or Focus. May be removed and interchanged with country-specific cords for travel. Contact Customer Service for country-specific cord part numbers.
	CD034-1	<b>DC Auto Adaptor Cable for FreeStyle 3 and Focus:</b> DC Adapter used to recharge and operate the FreeStyle 3 and Focus in an automobile.
	MI017-1	<b>FreeStyle &amp; Focus Carry-All Bag Strap:</b> Designed to allow the accessory bag to be carried over the shoulder. Attaches to MI320-1.
	MI284-1	<b>Backpack Harness:</b> Designed specifically for the FreeStyle 3 and FreeStyle 5 portable oxygen concentrators for easy transportation. MI284-1 includes solely the harness needed to convert either the FreeStyle or FreeStyle 5 carrying bag (MI406-1 or MI407-1) into a backpack.
	MI304-1	<b>FreeStyle Carrying Strap:</b> Allows the FreeStyle 3 or 5 to be carried over the shoulder. Can be attached to either MI406-1 or MI407-1.

# Accessories







Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

AirSep FreeStyle™		
Item Image	Item Number	Description
	MI320-1	<b>FreeStyle &amp; Focus Carry-All Bag:</b> Accessory bag designed specifically for easy transport of the FreeStyle and Focus portable oxygen concentrators and their accessories.
	MI406-1	<b>FreeStyle 3 Carrying Bag:</b> <i>NEW</i> Improved design. Case allows for easy transport of the FreeStyle by the carrying handle, or attachment of a shoulder strap or backpack harness. Bag also includes pockets to hold the External Power Cartridge for extended durations.
	<i>Country-Specific Power Cords Available:</i> PW022-1 (US) PW022-2 (Europe) PW022-3 (Australia) PW022-4 (UK)	<b>Universal Power Supply for FreeStyle 3:</b> A complete kit which includes everything needed to operate and recharge the FreeStyle 3 on either AC or DC power. Includes items CD034-1 and CD023-2.
	PW025-1	<b>Universal Power Supply (Brick Only) For FreeStyle:</b> Does not include any cords or cables.
AirSep FreeStyle™ 5		
Item Image	Item Number	Description
	20678526	<b>Wheelchair Straps:</b> Straps designed for easy attachment of concentrators to wheelchairs. Can be used with the Eclipse, eQuinox, FreeStyle, and FreeStyle 5.
	BT021-1S	<b>FreeStyle 5 Battery:</b> Internal battery for the FreeStyle 5 portable oxygen concentrator.
	CD023-2	<b>AC Power Cord:</b> Detachable hospital-grade US power supply that connects to the AC Power Supply of the LifeStyle, FreeStyle 3, FreeStyle 5, or Focus. May be removed and interchanged with country-specific cords for travel. Contact Customer Service for country-specific cord part numbers.







# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

AirSep FreeStyle™ 5		
Item Image	Item Number	Description
	MI284-1	<b>Backpack Harness:</b> Designed specifically for the FreeStyle 3 and FreeStyle 5 portable oxygen concentrators for easy transportation. MI284-1 includes solely the harness needed to convert either the FreeStyle or FreeStyle 5 carrying bag (MI406-1 or MI407-1) into a backpack.
	MI304-1	<b>FreeStyle Carrying Strap:</b> Allows the FreeStyle 3 or 5 to be carried over the shoulder. Can be attached to either MI406-1 or MI407-1.
	MI372-1	<b>FreeStyle 5 Carry-All Bag:</b> Accessory bag designed specifically for easy transport of the FreeStyle 5 portable oxygen concentrator and its accessories.
	MI407-1	<b>FreeStyle 5 Carrying Bag: NEW</b> Improved design. Case allows for easy transport of the FreeStyle 5 by the carrying handle, or attachment of a shoulder strap or backpack harness. Bag also includes pockets to hold the External Power Cartridge for extended durations.
	<i>Country-Specific Power Cords Available:</i> PW021-1 (US) PW021-2 (Europe) PW021-3 (Australia) PW021-4 (UK)	<b>AC Power Supply for FreeStyle 5:</b> Power Supply designed to recharge and operate the FreeStyle 5 in a standard wall outlet.
	PW027-1	<b>DC Auto Adapter and Power Supply for FreeStyle 5:</b> Power Supply designed to recharge and operate the FreeStyle 5 in an automobile.

# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

AirSep LifeStyle <sup>™</sup>		
Item Image	Item Number	Description
	BT007-1S	<b>LifeStyle Battery Pack:</b> Battery for LifeStyle.
	MI231-1	<b>LifeStyle Accessory Bag:</b> Accessory bag designed to hold up to six (6) LifeStyle batteries.
	PW005-1	<b>AC Power Supply for the LifeStyle:</b> Power Supply designed to recharge and operate the LifeStyle in a standard wall outlet.
	PW015-1	<b>DC Power Supply for the LifeStyle:</b> Power Supply designed to recharge and operate the LifeStyle in an automobile.

# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

SeQual Eclipse <sup>™</sup>		
Item Image	Item Number	Description
	20678526	<b>Wheelchair Straps:</b> Straps designed for easy attachment of concentrators to wheelchairs. Can be used with the Eclipse, eQuinox, FreeStyle, and FreeStyle 5.
	5010-SEQ	<b>Eclipse Travel Case:</b> Black, textured bag complete with four internal compartments for easy transportation of extra batteries, power supplies, etc. Easy attachment to the telescoping handle for those traveling with the Eclipse.
	5052-SEQ	<b>Eclipse Protective Cover:</b> Weather repellent cover for the Eclipse to aid in protection of the Eclipse from the elements during day-to-day use.
	5093-SEQ	<b>Travel Accessory Kit:</b> Perfect for Eclipse users that get up and go... and keep going! Includes the Eclipse Travel Case packed with an additional Power Cartridge and Desktop Charger.
	5220-SEQ	<b>Eclipse Wheelchair Pack:</b> Pack designed for easy attachment of the Eclipse to a wheelchair for easy transportation.
	5591-SEQ	<b>Universal Power Adapter:</b> Universal Power Adapter with adapters to connect into standard wall outlets in Europe, UK, New Zealand, North/South America, China, Japan, and Taiwan. Compatible with all portable and stationary oxygen concentrators.
	5941-SEQ	<b>Eclipse AC Power Supply:</b> US-style power supply designed to recharge the Eclipse in a standard wall outlet. Compatible with Eclipse 2, Eclipse 3, and Eclipse 5 models.

# Accessories

Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

SeQual Eclipse <sup>™</sup>		
Item Image	Item Number	Description
	5942-SEQ	<b>Eclipse DC Power Supply:</b> Power Supply designed to run the Eclipse in an automobile. Compatible with Eclipse 2, Eclipse 3, and Eclipse 5 models.
	5991-SEQ	<b>Eclipse Universal Cart with Telescoping Handle:</b> Pull-behind cart developed for easy transportation of the Eclipse transportable oxygen concentrator.
	7082-SEQ	<b>Eclipse Power Cartridge:</b> Compatible with all models of the Eclipse - Eclipse 1, Eclipse 2, Eclipse 3, and Eclipse 5.
	7104-SEQ	<b>Eclipse Accessory Bag:</b> Black, cloth bag designed for attachment to the telescoping handle for the Eclipse Universal Cart for easy transportation of personal items and travel items for day trips.
	7112-SEQ	<b>Eclipse Desktop Charger:</b> Tabletop charger designed specifically to recharge the Eclipse power cartridge. Includes charger and AC power supply.
SeQual eQuinox <sup>™</sup>		
Item Image	Item Number	Description
	20678526	<b>Wheelchair Straps:</b> Straps designed for easy attachment of concentrators to wheelchairs. Can be used with the Eclipse, eQuinox, FreeStyle, and FreeStyle 5.
	4823-SEQ	<b>eQuinox 12-Cell Battery Pack:</b> Standard battery pack for use with the SeQual eQuinox.

# Accessories








Concentrator Batteries, Power Supplies,  
Backpacks, Bags, Carts, Belts

SeQual eQuinox <sup>™</sup>		
Item Image	Item Number	Description
	4830-SEQ	<b>eQuinox DC Power Cable:</b> Cable to operate the eQuinox from an accessory outlet in a vehicle. Connects directly to the eQuinox.
	4836-SEQ	<b>eQuinox AC Power Supply (no cord):</b> Power supply used to operate the eQuinox from AC power.
	4880-SEQ	<b>eQuinox Cart:</b> Attaches to the eQuinox for easy transportation.
	4920-SEQ	<b>eQuinox Accessory Bag:</b> Bag designed to carry eQuinox batteries and power supplies along with personal items for day trips.
	4972-SEQ	<b>eQuinox 24-Cell Battery Pack:</b> Larger battery pack for extended use of the SeQual eQuinox.
	4980-SEQ	<b>eQuinox Desktop Charger:</b> Tabletop charger to recharge both the 12 and 24-cell battery packs from AC power. Includes charger and AC power supply.
	4997-SEQ	<b>eQuinox AC Power Cord:</b> US Power Cord that connects to the AC power supply (4836-SEQ) to operate and charge the eQuinox from a standard wall outlet.



# Accessories




## Concentrator Therapeutic Accessories

Item Image	Item Number	Description
	7116-SEQ	<b>SeQual Humidifier Adapter Kit:</b> Easy attachment allowing for the addition of a humidifier bottle to the SeQual eQuinox and Eclipse. <i>Please note: humidifier bottle is not included.</i>
	7568-SEQ	<b>SeQual Fingertip Oximeter:</b> Device designed to monitor oxygen saturation.
	CU002-1	<b>7 ft (2.1 m) Standard Nasal Cannula</b>
	CU002-4	<b>25 ft (7.6 m) Cannula/Tubing</b>
	CU004-3	<b>25 ft (7.6 m) Oxygen Tubing</b>
	F0025-1	<b>Oxygen Outlet Adaptor</b>
	HU003-1	<b>Humidifier Bottle:</b> For use with the Eclipse, eQuinox, NewLife Elite, and VisionAire concentrator models.




# Accessories

## Concentrator Testing and Service Equipment

### AirSep Focus<sup>™</sup>, FreeStyle<sup>™</sup>, and FreeStyle<sup>™</sup> 5






Item Image	Item Number	Description
	FI144-1	Air Intake Filter for FreeStyle
	FI194-1	Air Intake Filter for FreeStyle 5
	KI406-1	<b>AirSep Focus and FreeStyle 5 Hour Meter Kit:</b> USB Kit that allows you to check and reset the hours on the AirSep Focus and FreeStyle 5 via a computer.

### AirSep NewLife<sup>®</sup> and VisionAire<sup>™</sup> Families

Item Image	Item Number	Description
	FI002-1	Gross Particle Filter for NewLife Elite, Intensity, and Intensity 10
	KI036-1	<b>AirSep Test Pressure Gauge:</b> Gauge used to test the outlet and operating pressures of all AirSep Stationary Concentrators. Compatible with the tubing and fittings listed in this section.
	KI257-1	<b>AirSep NewLife Operating Pressure Test Fitting:</b> Fitting to allow connection of gauge (KI036-1) for testing the operating pressure of the NewLife Elite and Intensity.

# Accessories



## Concentrator Testing and Service Equipment

AirSep NewLife <sup>®</sup> and VisionAire <sup>™</sup> Families		
Item Image	Item Number	Description
	KI257-2	<b>AirSep VisionAire Operating Pressure Test Fitting:</b> Fitting to allow connection of gauge (KI036-1) to VisionAire for testing operating pressure.
	MI161-2 (NewLife Elite and NewLife Intensity 10) <i>For 120V NewLife Elite units, order MI161-1</i>	Optional Felt Filter for NewLife Concentrators
	MI161-3 (NewLife Intensity)	
	TU131-1	<b>AirSep NewLife Intensity 8 Outlet Pressure Tubing:</b> Tubing used to test the outlet pressure of the NewLife Intensity 8.
	TU131-2	<b>AirSep NewLife Intensity 10 Outlet Pressure Tubing:</b> Tubing used to test the outlet pressure of the NewLife Intensity 10.




# Accessories

## Concentrator Testing and Service Equipment

### SeQual Eclipse<sup>™</sup>

Item Image	Item Number	Description
	5022-SEQ	<b>PM Kit for the SeQual Eclipse 1, 2, 3, and 5:</b> Contains Air Intake Filter, Compressor Intake Filter, HEPA filter for either style product tank, 9V battery (with foam strip), and necessary zip ties.
	5535-SEQ	<b>SeQual EDAT:</b> Eclipse Data Acquisition Software. Real-Time diagnostic and quick-check computer program for all models of the SeQual Eclipse. Part number includes software, blue USB key, and necessary cable.

### SeQual eQuinox<sup>™</sup>

Item Image	Item Number	Description
	4770-SEQ	Oxygen Outlet Filter
	4845-SEQ	Air Intake Filter
	9759-SEQ	Compressor Inlet Filter for the eQuinox

# Liquid Oxygen Accessories



Discover your  
perfect **fit** from CAIRE.




# Accessories

## Liquid Oxygen Portable Backpacks

Companion <sup>™</sup>		
Item Image	Item Number	Description
	10018138	<b>C1000/C1000T Backpack:</b> Backpack is designed to fit Companion C1000 and C1000T portables. Contains detachable single strap for carrying on the shoulder as well as two backpack straps.
HELiOS <sup>™</sup>		
Item Image	Item Number	Description
	069209	<b>HELiOS Backpack H850:</b> Designed and sized specifically for the HELiOS Marathon (H850). Equipped with moisture pad to collect condensation. Can also be converted to a hip pack.
Hi Flow Stroller <sup>®</sup> , Spirit <sup>™</sup> , and Sprint <sup>™</sup>		
Item Image	Item Number	Description
	13214825	<b>Backpack:</b> Compatible with Spirit 600, Spirit 1200, Hi Flow, Stroller, and Sprint portables.

# Accessories





## Liquid Oxygen Portable Hip Packs

<b>Companion<sup>™</sup></b>		
Item Image	Item Number	Description
	B-778625-00	<b>C500/C550 Hip Pack:</b> Pack for the Companion C500 and C550 portables designed to be worn on the hip.
<b>HELiOS<sup>™</sup></b>		
Item Image	Item Number	Description
	14951835	<b>Large Belt for H300 Hip Pack:</b> Large replacement belt for the H300 hip pack only (B-701654-00). 48-81 inches (121.9-205.7 cm)
	B-701654-00	<b>H300 Hip Pack:</b> Pack for the HELiOS Plus (H300) designed to be worn on the hip.

# Accessories


## Liquid Oxygen Portable Hip Packs

### Spirit<sup>™</sup>

Item Image	Item Number	Description
	11843361	<b>Spirit 300 Belt Pack:</b> Pack for the Spirit 300 that is designed to be worn on the hip.
 Attaches to existing straps here	11843370	<b>Belt Pack Side Pouch:</b> Accessory pouch to be used with the Spirit 300 belt pack (#11843361). Slides into the waist strap of the belt pack to allow for additional storage of accessories. 7 in x 4 in (17.8 cm x 20.3 cm) with a zipper top.
	11843388	<b>Belt Pack Side Pouch:</b> Accessory pouch to be used with the Spirit 300 belt pack (#11843361). Slides onto the waist strap of the belt pack to allow for additional storage of accessories. 7 in x 8 in (17.8 cm x 20.3 cm) with a zipper top.
	11843396	<b>Spirit 300 Belt Pack Waist Extension:</b> Adjustable belt extension for the Spirit 300 belt pack. Compatible with part number #11843361 only.





# Accessories

## Liquid Oxygen Portable Universal Cart

Item Image	Item Number	Description
	20773616	<b>Universal Liquid Oxygen Portable Cart:</b> Two-wheeled pull behind cart with a Velcro strap to secure portables. Compatible with Stroller, Hi-Flow, HELiOS Marathon (H850), and Companion 1000/1000T portables.

# Accessories



## Liquid Oxygen Reservoir Rollerbase Assemblies

Item Image	Item Number	Description
	CE Marked: 10855678	<b>Companion Reservoir Rollerbase:</b> 14" roller base designed to easily roll and transport the Companion C21, C31, and C41 reservoirs.
	Non CE Marked: 10585866	
	CE Marked: 14880387	<b>HELIOS Reservoir Rollerbase:</b> 15" roller base designed to easily roll and transport the HELiOS H36, H46, U36, and U46 reservoirs.
	Non CE Marked: 14880803	
	CE Marked: 11682262	<b>Liberator 10 Rollerbase:</b> 12" roller base designed to easily roll and transport the Lib 10 reservoir.
	Non CE Marked: 10585858	
	CE Marked: 10855678	<b>Liberator 20, 30, 37, 45 Rollerbase:</b> 15" roller base designed to easily roll and transport the Liberator 20, 30, 37, and 45 reservoirs.
	Non CE Marked: 10585866	
	CE Marked: 11363069	<b>Liberator 60 Rollerbase:</b> 16" roller base designed to easily roll and transport the Liberator 60 reservoirs.
	Non CE Marked: 11194074	
<b>NOTE: CE marked rollerbases are approved for sell and distribution in Europe. They contain 2 wheels that contain brakes, and 3 wheels that do not. Non CE marked reservoirs contain 5 wheels without brakes.</b>		
	14880379	<b>CE Marked Casters:</b> Replacement set of 3" casters for CE marked reservoir roller base assemblies. Include 5 casters total. 2 with brakes and 3 without brakes.
	14880361	<b>Non CE Marked Casters:</b> Replacement set of 3" casters for non CE marked reservoir roller base assemblies. Includes 5 casters without brakes.








# Accessories

## Nasal Cannulas

Item Image	Item Number	Description
	4 ft (1.2 m): B-701931-00	<b>Dual Lumen Cannula:</b> Demand nasal cannula for use with HELiOS portables. Sense and delivery in separate nostrils (oxygen will only be delivered through 1 nostril).
	7 ft (2.1 m): B-701930-00	
	5 ft (1.5 m): 6-778058-00	<b>Dual Lumen Cannula:</b> Demand nasal cannula for use with HELiOS portables. Sense and delivery in both nostrils with a divided tube.
	7 ft (2.1 m): 6-778057-00	
	5 ft (1.5 m): 10035468	<b>Dual Lumen Cannula:</b> Demand nasal cannula for use with HELiOS portables. Sense and delivery in both nostrils using two tubes.
	7 ft (2.1 m): 10035467	
	4 ft (1.2 m): 20574985	<b>Dual Lumen Cannula:</b> Demand nasal cannula for use with HELiOS portables. Sense and delivery in both nostrils with a single tube.
	5 ft (1.5 m): 20574977	
	7 ft (2.1 m): 2057497	
	5408-SEQ	<b>Standard Nasal Cannula:</b> Single lumen 7 ft (2.1 m) cannula for continuous flow portables and reservoirs.






# Accessories

## Transfilling Accessories

Item Image	Item Number	Description
	3910486	<b>Dewar Cap:</b> Used to cover the opening of Liberator and Companion reservoir dewars when the manifold has been removed.
 	4 ft (1.2 m): 9713139	<b>Transfer Hose:</b> Standard oxygen service transfer line with swivel nuts used for transfilling liquid oxygen reservoirs. One end connects to the transfill adaptor and the other end connects to the oxygen source.
	6 ft (1.8 m): 9713119	
	8 ft (2.4 m): 10546550	
	10 ft (3.0 m): 10565161	
	12 ft (3.7 m): 10562411	
	10678157	<b>CAIRE Top Fill Transfer Line Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all top fill reservoirs with blue male PB fittings.
	10768831	<b>Cryo-2 Transfer Line Adaptor:</b> Connects to the end of a transfer line. Used for transfilling all top fill reservoirs with a Cryo-2 fitting.

# Accessories

## Transfilling Accessories

Item Image	Item Number	Description
	10769294	<b>Penox Female Transfer Line Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all top fill (TF) reservoirs with a Penox top fill fitting.
	10897958	<b>Dual Fill Transfill Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all top fill (TF) reservoirs with a blue male PB fill connector and side fill (SF) Liberator reservoirs.
	10917963	<b>Linde TF &amp; CAIRE SF Dual Fill Transfill Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into top fill (TF) reservoirs with Linde fittings and side fill (SF) Liberator units.
	10917971	<b>CAIRE SF &amp; Pulseair TF Dual Fill Transfill Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all side fill (SF) Liberator units and all top fill reservoirs with Pulseair fittings.
	10917980	<b>CAIRE SF &amp; Penox TF Dual Fill Transfill Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all side fill (SF) Liberator units and all top fill (TF) reservoirs with Penox fittings.






# Accessories

## Transfilling Accessories

Item Image	Item Number	Description
	10939601	<b>CAIRE TF &amp; Linde TF Dual Fill Transfill Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all top fill (TF) reservoirs with a blue male PB fill connector and top fill reservoirs with Linde fittings.
	97202005	<b>Reservoir Vent Valve Wrench:</b> Vent valve wrench for all Liberator reservoirs manufactured before July 2012. This vent valve wrench may also be used on all CAIRE reservoirs manufactured after July 2012.
	97212023	<b>Male Side Fill Transfer Line Adaptor:</b> Connects to the end of a transfer line. Used for transfilling into all Side Fill (SF) Liberator reservoirs.
	97217007	Fitting Adapter 90 D
	97404564	<b>Swivel Connector:</b> Male end connects to transfer hose and female end connects to transfill adaptor. This allows the transfill adaptor to rotate freely.

# Accessories

## Transfilling Accessories

Item Image	Item Number	Description
	97406555	<b>Super Flex Hose:</b> 6 ft (1.8 m) oxygen service transfer line without swivel nuts. Used for filling liquid oxygen reservoirs. One end connects to the transfill adaptor and the other end connects to the oxygen source.
	B-701732-00	<b>Pressure Gauge Assembly:</b> Attaches to the flow control valve (FCV) outlet of a reservoir during transfill. Used to monitor and maintain the internal pressure of the reservoir during a fill.
	B-775182-00	<b>Reservoir Vent Valve Wrench:</b> Vent valve wrench for all HELIOS and Companion reservoirs manufactured before July 2012. This vent valve wrench may also be used on all CAIRE reservoirs manufactured after July 2012.
	CA400004	<b>Inline Filter:</b> Replacement filter for all transfill line adaptors.
	CA405311	<b>Stroller Vent Key:</b> Vent valve key for G2 and G2.5 Stroller portable units (older tan hard cases).








# Accessories

## Testing Accessories

Item Image	Item Number	Description
	10679862	<b>Female Top Fill Pneumatic Test Adaptor:</b> Connects to a gaseous oxygen or nitrogen source (can use pneumatic hose #97405279) for gaseous testing of reservoirs with blue male top fill connectors.
	10679871	<b>Male Top Fill Pneumatic Test Adaptor:</b> Connects to a gaseous oxygen or nitrogen source (can use pneumatic hose #97405279) for gaseous testing of portable units with blue female top fill connectors.
	10995620	<b>Erie Liter Meter, 6-15 LPM:</b> Used to measure continuous flow from reservoirs and portables between 6-15 LPM.
	11781228	<b>Pulsemeter:</b> Used to test the accuracy of pulse dosing with each breath of HELiOS and Spirit portables.
	13095635	<b>Solder 60/40 Lead Free:</b> Used to solder new wire harness onto Liberator reservoirs and Spirit/Sprint/Stroller portables with a G4 contents indicator meter.
	13188170	QDV Assembly for Gas Analyzer







# Accessories

## Testing Accessories

Item Image	Item Number	Description
	13329091	<b>G4 Capacitance Kit:</b> Used with the Capacitance Meter to test the capacitance of the G4 meter of the Liberator reservoir or Sprint/Stroller portable.
	97112026	<b>Female Side Fill Pneumatic Test Adaptor:</b> Connects to a gaseous oxygen or nitrogen source (can use pneumatic hose #97405279) for gaseous testing of portable units with side fill (SF) connections.
	97200076	<b>Erie Liter Meter, 0-8 LPM:</b> Used to measure continuous flow from reservoirs and portables between 0-8 LPM.
	97212021	<b>Male Side Fill Pneumatic Test Adaptor:</b> Connects to a gaseous oxygen or nitrogen source (can use pneumatic hose #97405279) for gaseous testing of reservoir units with side fill (SF) connections.
	97403015	<b>Capacitance Meter:</b> Meter used to measure capacitance of the manifold on Liberator reservoirs and Sprint/Stroller portables with meters. Important for troubleshooting issues with the G4 content indicator.






# Accessories

## Testing Accessories

Item Image	Item Number	Description
	97403577	<b>0-60 psi (0-4.14 bar) Pressure Gauge:</b> Used for testing or verifying internal pressure of LOX portables and reservoirs.
	97405279	<b>Pneumatic Hose:</b> Hose has DISS fittings on both ends. Used to connect a gaseous source of oxygen or nitrogen to test adaptors for gaseous testing of reservoirs and portables.
	97405590	<b>Lip Seal Tool:</b> Tool used to remove the locking mechanism on Side Fill Liberator reservoirs so that the lip seal can be removed.
	AAA0645	<b>Soldering Flux:</b> Used to clean the capacitance probe when soldering a new wire harness onto the probe. Use on Liberator reservoirs and Sprint/Stroller/Hi Flow portables with a G4 meter.
	B-701731-SV	<b>Reservoir Pressurizing Fixture:</b> Used with all top fill reservoirs to pneumatically pressurize the unit. Used for testing purposes and also for emptying and purging moisture from a reservoir.
	B-775239-00	<b>Krytox Lubricant:</b> Lubrication used on fittings to help prevent cross-threading. Used on HELIOS and Companion reservoirs and portables.







# Accessories

## Testing Accessories

Item Image	Item Number	Description
	B-775272-00	<b>Liquid Snoop (8 oz):</b> Used for detecting leaks on all portables and reservoirs.
	B-778202-00	<b>Portable Test Fixture:</b> Used to hold portable units during function testing. Prevents the portables from lying flat and makes easy access to all components. Especially useful for servicing HELiOS portables with the case removed.
 	B-778210-00	<b>Jet Venturi Assembly:</b> Essential component for testing the conservers of the H300 and H850. One end connects to a gaseous source and the other end connects to the bottom barb of the HELiOS portable to trigger demand flow.
	B-778894-00	<b>Liquid Snoop (1 gal):</b> Used for detecting leaks on all portables and reservoirs.




# Accessories

## Testing Accessories

Item Image	Item Number	Description
	CA110104	<b>Side Fill Lip Seal:</b> Plastic fitting that seals the mating male and female QDV's during the filling process from reservoir to portable.
	CA200071	<b>Christo-Lube:</b> Lubrication used on fittings to help prevent cross-threading. Used on Liberator, Stroller, Sprint, and Spirit product lines.
	CA404838	O-Ring For Male 1-Pin ODV SF
	CA406310	<b>TEFLON Tape:</b> Tape used on threaded connections to prevent leaks at the connections. Used on all reservoirs and portables.
	CA406354	Swivel Connector Kit
	CA406506	<b>Static Seal:</b> Plastic fitting used on top fill QDV's that seals the mating male and female QDV's during the filling process from reservoir to portable.

# Accessories

## General Accessories

Item Image	Item Number	Description
	B-701655-FLO	<p><b>HELiOS Reservoir External FCV:</b> 0-10 LPM Flow Control Valve for patients who wish to breathe directly from the HELiOS reservoir.</p> <p><i>NOTE: This item is not compatible with EU models manufactured before December 2008. These units use #10000554 as an external FCV.</i></p>
	B-701656-SV	<p><b>HELiOS Oxygen Supply Line:</b> 50 ft (15.2 m) line used to connect the HELiOS H300 or H850 to a reservoir. When connected, the patient can breathe the oxygen from the reservoir using the portable unit as a conserver.</p>
	B-776945-00	<p><b>DISS Adaptor:</b> Connects to the O2 outlet of a reservoir so that a patient can connect tubing or a cannula for direct breathing.</p>



# Accessories

## Flow Upgrades and Conversions

Item Image	Item Number	Description
	13440224	<b>Liberator 15LPM Upgrade Kit:</b> Designed to upgrade older 6 LPM and 10 LPM Liberator reservoirs to 0-15 LPM flow capabilities. Kit includes: Replacement warming coil, 0-15 LPM FCV, cable ties (x3), nuts (x2), and ferrules (x2)
 FCV Extension  Drain Tube  Elbow and Bracket  Condensate Bottle	13660403	<b>Liberator Accessory Kit:</b> Contains condensate bottle, condensate elbow and bracket, and the G4 flow control valve extension.
	Pediatric Decal: 14700411 Pediatric FCV: 14700235T	<b>Sprint/Stroller Pediatric Conversion:</b> Parts required to replace the flow control valve on G3 Liberator reservoir with a pediatric FCV capable of achieving these low rates: 0.05, 0.10, 0.20, 0.30, 0.40, 0.50, 0.75, 1.00, 1.25, 2.00 LPM.



# Accessories

## Flow Upgrades and Conversions

Item Image	Item Number	Description
 Decal 1470038	Pediatric Decal: 14700438	<b>Gen 3 Liberator Pediatric Conversion:</b> Parts required to replace the flow control valve on G3 Liberator reservoir with a pediatric FCV capable of achieving these flow rates: 0.05, 0.10, 0.20, 0.30, 0.40, 0.50, 0.75, 1.0, 1.25, 2.00 LPM
 FCV 14700243T	Pediatric FCV: 14700243T	
 Decal 1470038	Pediatric Decal: 14700438	<b>Gen 4 Liberator Pediatric Conversion:</b> Parts required to replace the flow control valve on G3 Liberator reservoir with a pediatric FCV capable of achieving these flow rates: 0.05, 0.10, 0.20, 0.30, 0.40, 0.50, 0.75, 1.0, 1.25, 2.00 LPM
 FCV 14700243T	Pediatric FCV: 14700243T	
 CA405938	FCV Extension: CA405938	
 CA406438	Locking Plate: CA406438	
 10495289T	SRV: 10495289T	






# Accessories

## Flow Upgrades and Conversions

Item Image	Item Number	Description
	B-701362-SV	<b>Companion Reservoir Drain Bottle Kit:</b> Contains the drain bottle, drain tube, and bracket to fit all Companion reservoirs. The drain bottle attaches to the side of the reservoir and is used to collect condensation after frost on the coils as it melts.
	B-701432-00	<b>Supply Line Tubing Extension:</b> 50 ft (15.2 m) extension tubing for the HELiOS oxygen supply line (B-701656-SV).

# Accessories

## CAIRE Generation 4 Meter Upgrades

Item Image	Item Number	Description
	13199653	<b>G2 to G4 Upgrade Kit for Portables and Reservoirs:</b> Kit to upgrade older Generation 2 Stroller portables and Liberator reservoirs to the current Gen 4 contents indicator meter. Kit contains: G4 meter, wire harness adaptor, and replacement Stroller top cover. The Stroller top cover will not be used if upgrading a reservoir.
	13199653	<b>G4 High Flow Upgrade Kit:</b> Kit to upgrade older Hi Flow portables (tan hard case) to the current Gen 4 contents indicator meter. Kit contains: G4 meter, wire harness adaptor, and replacement top cover.
	13199661	<b>G3 to G4 Stroller/Sprint Upgrade Kit:</b> Kit to upgrade older Generation 3 Stroller and Sprint portable units to the current Gen 4 contents indicator meter. Kit contains: G4 meter and wire harness adaptor. <i>Must order #12932526 (replacement dome cover) with the kit to complete the upgrade.</i>
	13219221	<b>G3 to G4 Liberator Upgrade Kit:</b> Kit to upgrade older Generation 3 Liberator reservoirs to the current Gen 4 contents indicator meter. Kit contains: G4 meter, ground wire, and wire harness adaptor.
	13579584	<b>G3 to G4 Low Loss/EMS 20 Upgrade Kit:</b> Kit to upgrade older Low Loss and EMS 20 reservoirs to the current Gen 4 contents indicator meter. Kit contains: G4 meter, wire harness adaptor, ground cable, and telemetry bracket. Must be ordered with the harness adaptor nut (#13717262).
<b>NOTE: All G4 meters received as part of an Upgrade Kit are not calibrated. These will need to be calibrated per the included instructions prior to use.</b>		

# Accessories

## Spirit Upgrade Kits

Item Image	Item Number	Description
	12937052	<b>Simulated Leather Retro Kit 300 Side Fill</b>
	12937061	<b>Simulated Leather Retro Kit 300 Top Fill</b>
	13029030	<b>Battery Holder Upgrade</b>
	13106251	<b>Knob Adaptor Repair Kit</b>

All illustrations and specifications contained herein are based on the latest product information at the time of publication. CAIRE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

# Accessories Catalog



Discover your  
perfect**fit** from CAIRE.

CAIRE, Inc.  
2200 Airport Industrial Drive, Suite 500  
Ball Ground, GA 30107  
Ph 770-721-7759 • Toll-Free 1-800-482-2473 Fax 888-932-2473

Chart BioMedical, Ltd.  
Unit 6, Ashville Way  
Wokingham RG41 2PL, United Kingdom  
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245